Laboratory Heating Mantles Instruction Manual

	<p>Read and understand this manual before operating or servicing this Laboratory Heating Mantle. Failure to understand how to safely operate this heater could result in an accident causing serious injury or death. Only qualified personnel should operate or service this heater.</p>
---	---

Language	Page
English	1
Spanish (Español).....	9
French (Français).....	17
German (Deutsch).....	25
Italian (Italiano).....	33

TABLE OF CONTENTS

Introduction.....	2
Safety Alert Symbol	2
Important Safety Instructions.....	3
Specifications	4
Prior To Installation.....	5
Installation	6
Troubleshooting Guide	8
Warranty Information	8

SAFETY ALERT SYMBOL

INTRODUCTION

Thank you for purchasing a BriskHeat[®] Laboratory Heating Mantle. Your heater is designed to provide a long and efficient service life with function, reliability, and safety in mind. For additional information concerning this, or other BriskHeat[®] products, please contact your local BriskHeat[®] distributor or contact us toll free (U.S. / Canada) at 1-800-848-7673 or 614-294-3376 (Worldwide).

The symbol above is used to call your attention to instructions concerning your personal safety. It points out important safety precautions. It means “**ATTENTION! Become Alert! Your Personal Safety is involved!**” Read the message that follows and be alert to the possibility of personal injury or death.

Immediate hazards which **WILL** result in severe personal injury or death.

Hazards or unsafe practices that **COULD** result in severe personal injury or death.

Hazards or unsafe practices that **COULD** result in minor personal injury or property damage.

SAVE THESE INSTRUCTIONS!

Additional copies of this manual are available upon request.

IMPORTANT SAFETY INSTRUCTIONS

⚠ DANGER

A person who has not read and understood all operating Instructions is not qualified to operate this product.

⚠ DANGER

- Do not immerse or spray any component of the control system with liquid.
- Keep volatile or combustible material away from controller when in use.
- Use heater only in approved locations.
- Keep sharp metal objects away from heater.

Failure to observe these warnings may result in electric shock, risk of fire, and personal injury.

⚠ WARNING

End-User Must Comply to the Following:

- Only qualified personnel are allowed to connect electrical wiring.
- Disconnect all supply power at the source before making any power connections.
- All electrical wiring must follow local electrical codes and highly recommend following NEC Article 427.
- The person who performs the final installation / wiring must be qualified for this work.
- The end-user is responsible for providing a suitable disconnecting device.
- The end-user is responsible for providing suitable electrical protection device. It is highly recommended that a ground fault circuit breaker be used.

Failure to observe these warnings may result in personal injury or damage to the heater.

Agency Approvals

Approvals valid only when installed in accordance with all applicable instructions, codes, and regulations.

⚠ CAUTION

- Never handle the heater while it is in operation; always disconnect the heater from the power source and allow to cool prior to handling.
- Inspect heater before use.
- Never operate heater without a temperature control device.
- Do not wrap the heater over itself.
- If spillage of foreign matter onto heater occurs, disconnect from power source and clean after heater is allowed to cool.
- Never operate a heater without an appropriate heat sink (device being heated is considered a heat sink).
- Do not operate heater above rated temperature value.
- Fasten heater to device using approved methods only.
- Do not repair damaged or faulty heaters.
- Do not crush or apply severe physical stress on heater or cord assembly.
- Disconnect heater when not in use.

Failure to observe these warnings may result in personal injury or damage to the heater.

SPECIFICATIONS

		HM Lower Hemispherical	HM Upper Hemispherical	HM Spherical	HM Table Top
Product Highlights		Ideal for use with basket ring stand.	Opening holds 1, 2 or 3 neck flasks on ring stand. Ideal for preventing vapor condensation in upper half of flask.	Opening holds 1, 2 or 3 neck flasks on ring stand, Zipper fastens upper and lower halves together.	Self-supporting can be placed directly on counter top. Easy viewing of flask contents. Stackable for easy storage.
Maximum Exposure Temperature		450°C (842°F)	450°C (842°F)	450°C (842°F)	482°C (900°F)
Part No. 120 VAC	50 ml 100 ml 125 ml 200 ml 250 ml 500 ml 1000 ml 2000 ml 3000 ml 5000 ml 12000 ml	HM0050-HS1 HM0100-HS1 HM0125-HS1 HM0200-HS1 HM0250-HS1 HM0500-HS1 HM1000-HS1 HM2000-HS1 HM3000-HS1 HM5000-HS1 HM12000-HS1	HM0250-TS1 HM0500-TS1 HM1000-TS1 HM2000-TS1 HM3000-TS1 HM5000-TS1	HM0500-SS1 HM1000-SS1 HM2000-SS1	HM0050-VF1 HM0100-VF1 HM0250-VF1 HM0500-VF1 HM1000-VF1 HM2000-VF1 HM3000-VF1 HM5000-VF1
Part No. 230 VAC	50 ml 100 ml 125 ml 200 ml 250 ml 500 ml 1000 ml 2000 ml 3000 ml 5000 ml 12000 ml	HM0050-HS2 HM0100-HS2 HM0125-HS2 HM0200-HS2 HM0250-HS2 HM0500-HS2 HM1000-HS2 HM2000-HS2 HM3000-HS2 HM5000-HS2 HM12000-HS2	HM0250-TS2 HM0500-TS2 HM1000-TS2 HM2000-TS2 HM3000-TS2 HM5000-TS2	HM0500-SS2 HM1000-SS2 HM2000-SS2	HM0050-VF2 HM0100-VF2 HM0250-VF2 HM0500-VF2 HM1000-VF2 HM2000-VF2 HM3000-VF2 HM5000-VF2

⚠ IMPORTANT: Temperature controller is required for this product.

PRIOR TO INSTALLATION

1. Check for suspected damage to the heater like rips, punctures, etc.
2. Confirm voltage / wattage rating of heater is appropriate for power supply device or temperature controller. The heater is designed to work at a specific voltage to heat properly. Failure to operate the heater at the specified rating could result in overheating of the heating blanket, the surface being heated, or personal injury. Also note lower-than-rated voltage will result in under-heating.

To Avoid Electrical Shock, Always:

- Use an electrical outlet that operates with a fuse or a circuit breaker and a ground fault circuit interrupt (GFCI).
- Disconnect the mantle from the power supply prior to maintenance and servicing.

To Avoid Personal Injury:

- Do not use in the presence of flammable or combustible materials, fire or explosion may result. This device contains components which may ignite such materials.
- Do not operate in damp or wet locations.
- Refer servicing to qualified personnel only.
- Flasks and vessels should be handled with care. Appropriate safety clothing, glasses, gloves, and coats should be worn when operating all heating mantles.

To Avoid Damage:


- Never plug directly into a power source. Your heating mantle MUST be operated with a temperature controller and with an appropriate size flask, filled with fluid.
- Never operate heating mantle without fluid within the flask.
- Fluid should not be allowed to come in contact with the cavity of the heating mantle.
- Never operate the mantle at temperatures above maximum temperature exposure.

INSTALLATION

Your Heating Mantle is a Fiberglass-insulated heater designed for heating laboratory vessels. The mantles are specifically designed for glass flasks, however, with appropriate temperature control and product type, metal or plastic flasks can be used. Only use round bottom flasks of the size corresponding to the mantle size.

Spherical and hemispherical mantles should be used with properly sized supports. The manufacturer does not assume responsibility for mantles damaged as a result of inadequate support.

1. Install heater around properly sized flasks.
2. For proper temperature control, secure the thermocouple, RTD, or other sensing device, vertically along the lead wire seam of the blanket with one strip of adhesive tape so that it is sandwiched firmly in-between the heater and glass flask.
 - Sensor must be in close proximity to the heater to prevent overheating.
 - Sensor should be connected in accordance with instructions provided with temperature controller.
3. Ensure heater is in contact with the surface to be heated. The heater must fit tightly around the glass flask. Minimal air gaps between the heater and the flasks provide the best heat transfer and best heater life.
4. Once the mantle has been properly setup:
 - 120VAC: Plug the heater into an acceptable temperature controller. Use a variable transformer, manual control, or automatic controller to prevent overheating.
 - 230VAC: The two wires (same color) must be connected to the temperature control device directly or to a connector that is mated to a receptacle connected to the temperature control source. These heaters are non-polarized, allowing either wire to be connected to neutral or hot.

All electrical connections must be made by qualified personnel and in accordance with all applicable codes and regulations.

Special Instructions for Spherical or Hemispherical Mantles with multiple circuits

The 1- and 2-liter spherical mantles have two circuits – one in the lower half and one in the upper half. Each circuit should be controlled with a suitable control device. The lower half's circuit furnishes heat for boiling the liquid contents, while the upper half's circuit prevents condensation of the vapors. It may be necessary to operate the upper half's circuit at no more than 50% of full power. For low boiling liquids, the upper half's circuit need not be used. This will prevent superheating the vapors, and overheating of the mantle's circuits.

The 12-liter hemispherical (bottom half) mantle also has two circuits. When a flask is more than half full of liquids both circuits may be operated at the rated voltage. When the liquid level falls below the halfway mark in the flask, power to the upper circuit of the mantle should be reduced to 50% or less. This will prevent superheating of the vapors, and overheating of the mantle's circuits.

Bake-out Procedure for New Heating Mantles

On the initial heat-up of the mantles, a slight odor may be detected and some discoloration will occur in the cavity area. The discoloration of the cavity area has no effect on the operation of the mantle. If the heated vessel becomes discolored, it can be cleaned with water. To bake-out your mantle prior to use, connect the mantle to your temperature controller and fill a flask corresponding to the size of the mantle half full of water. Place the flask in the mantle cavity and allow the water to come to a boil and continue boiling for approximately 15 minutes. The mantle can be operated at full rated voltage during bake-out.

Maintenance and Servicing

Your Heating Mantle is designed to provide a long service life. Overheating, contamination, and misuse will greatly reduce the life of the mantle.

If the flask or vessel breaks or fluid spills into the cavity, immediately disconnect the mantle from the power source.

TROUBLESHOOTING GUIDE

Please read this guide prior to contacting BriskHeat[®]. This guide is designed to answer the most commonly asked questions. If you are unable to identify the problem or need additional assistance, please contact your local distributor/ representative or us at **1-800-848-7673**, **614-294-3376**, or **bhtsales1@briskheat.com**.

PROBLEM	POSSIBLE CAUSE	CORRECTIVE ACTION
Heating Mantle does not heat	Not connected to power supply. No Power in circuit. Defective Control. Element burned out.	Check Mantle connection. Check power supply circuit. Repair or replace Control. Replace Mantle.
Heating Mantle outgasses	Has not been baked-out. Spillage of fluid.	Refer to Bake-Out Procedure for new mantles. Disconnect from power supply and allow mantle to dry out completely before use.

WARRANTY INFORMATION

BriskHeat warrants to the original purchaser of this product for the period of eighteen (18) months from date of shipment or twelve (12) months from date of installation, whichever comes first. BriskHeat's obligation and the exclusive remedy under this warranty shall be limited to the repair or replacement, at BriskHeat's option, of any parts of the product which may prove defective under prescribed use and service following BriskHeat's examination, is determined by BriskHeat to be defective. The complete details of the warranty can be found online at www.briskheat.com or by contacting us at 1-800-848-7673 (toll free, U.S. / Canada) or 1-614-294-3376 (Worldwide).

4800 Hilton Corporate Dr Columbus, OH 43232
Toll Free: 800-848-7673
Phone: 614-294-3376
Fax: 614-294-3807
Email: bhtsales1@briskheat.com

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Corporation